

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS
FORMULATED STANDARDS SETS**

Subcategory: CHLORINE AND SULFUR STANDARDS

Table: CHLORINE AND SULFUR STANDARDS

| Product Code | Product | Description | Price \$ Set | Individual (100 ml) | Individual (1 L) |
|---------------------|-------------------------------------|--|---------------------|----------------------------|-------------------------|
| IMCLSMO10 | CHLORINE AND SULFUR IN MINERAL OIL. | Chlorinated hydrocarbon and Di-N-butyl sulfide in mineral oil. Conc. Range: Cl:0-1.00 wt% S:0-1.00 wt% 10 standards per set, 100mls each | \$815.00 | \$110.00 | \$575.00 |
| IMCLSWO10 | CHLORINE AND SULFUR IN WASTE OIL. | Chlorinated hydrocarbon and Di-N-butyl sulfide in waste oil. Conc. Range: Cl:0-1.00 wt% S:0-2.50 wt% 10 standards per set, 100mls each | \$835.00 | \$125.00 | \$595.00 |

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS
FORMULATED STANDARDS SETS**

Subcategory: CHLORINE AND SULFUR STANDARDS

Table: IMCLSMO10

| Standard No. | Sulfur Wt % | Chlorine Wt. % |
|---------------------|------------------------|---------------------------|
| 1 | 0.0000 | 0.0000 |
| 2 | 0.1000 | 0.8000 |
| 3 | 0.6000 | 0.5000 |
| 4 | 0.2000 | 0.3000 |
| 5 | 0.7000 | 0.2000 |
| 6 | 0.3000 | 0.6000 |
| 7 | 0.4000 | 0.1000 |
| 8 | 0.5000 | 0.4000 |
| 9 | 0.8000 | 1.0000 |
| 10 | 1.0000 | 0.0000 |

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS
FORMULATED STANDARDS SETS**

Subcategory: CHLORINE AND SULFUR STANDARDS

Table: IMCLSWO10

| Standard No. | Sulfur Wt % | Chlorine Wt. % |
|---------------------|------------------------|---------------------------|
| 1 | 0.5000 | 1.0000 |
| 2 | 2.5000 | 0.0200 |
| 3 | 0.5000 | 0.6500 |
| 4 | 1.0000 | 0.2000 |
| 5 | 1.5000 | 0.5000 |
| 6 | 1.7500 | 0.8000 |
| 7 | 1.2500 | 0.0500 |
| 8 | 2.0000 | 1.0000 |
| 9 | 2.2500 | 0.3800 |
| 10 | 0.7500 | 0.1000 |