

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: LEAD STANDARDS**

**Table: LEAD STANDARDS**

<b>Product Code</b>	<b>Product</b>	<b>Method</b>	<b>Description</b>	<b>Price \$ Set</b>	<b>Individual (100 ml)</b>
IMBIIS	BISMUTH INTERNAL STANDARD	ASTM 5059	Internal Standard Conc: 0.793g Bi/L	\$350.00	\$90.00
IMPBG12	LEAD IN GASOLINE <b>FLAMMABLE UN1203</b>	ASTM D5059, PARTS A & C	This set is made up of PBG7A and PBG7C 12 standards per set, 100mls each	\$975.00	\$85.00
IMPBG7A	LEAD IN GASOLINE <b>FLAMMABLE UN1203</b>	ASTM D5059, PART A	7 standards per set, 100mls each	\$515.00	\$85.00
IMPBG7C	LEAD IN GASOLINE <b>FLAMMABLE UN1203</b>	ASTM D5059, PART C	7 standards per set, 100mls each	\$515.00	\$85.00
IMPBI07A	LEAD IN ISOOCTANE <b>FLAMMABLE UN1262</b>	ASTM D5059, PART A	7 standards per set, 100mls each	\$515.00	\$85.00
IMPBI012	LEAD IN ISOOCTANE <b>FLAMMABLE UN1262</b>	ASTM D5059, PARTS A & C	This set is made up of PBIO7A and PBIO7C 12 standards per set, 100mls each	\$975.00	\$85.00
IMPBI07C	LEAD IN ISOOCTANE <b>FLAMMABLE UN1262</b>	ASTM D5059, PART C	7 standards per set, 100mls each	\$515.00	\$85.00
IMPBT12	LEAD IN TOLUENE <b>FLAMMABLE UN1294</b>	ASTM D5059, PARTS A & C	This set is made up of PBT7A and PBT7C 12 standards per set, 100mls each	\$975.00	\$85.00
IMPBT7A	LEAD IN TOLUENE <b>FLAMMABLE UN1294</b>	ASTM D5059, PART A	7 standards per set, 100mls each	\$515.00	\$85.00
IMPBT7C	LEAD IN TOLUENE <b>FLAMMABLE UN1294</b>	ASTM D5059, PART C	7 standards per set, 100mls each	\$515.00	\$85.00

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: LEAD STANDARDS**

**Table: IMPBG7A**

<b>Standard No.</b>	<b>Lead (Grams/Gal.)</b>
1	0.000
2	0.100
3	1.000
4	2.000
5	3.000
6	4.000
7	5.000

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: LEAD STANDARDS**

**Table: IMPBG7C**

<b>Standard No.</b>	<b>Lead (Grams/Gal.)</b>
1	0.000
2	0.001
3	0.005
4	0.010
5	0.050
6	0.100
7	0.300

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: LEAD STANDARDS**

**Table: IMPBI07A**

<b>Standard No.</b>	<b>Lead (Grams/Gal.)</b>
1	0.000
2	0.100
3	1.000
4	2.000
5	3.000
6	4.000
7	5.000

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: LEAD STANDARDS**

**Table: IMPBIO7C**

<b>Standard No.</b>	<b>Lead (Grams/Gal.)</b>
1	0.000
2	0.001
3	0.005
4	0.010
5	0.050
6	0.100
7	0.300

<b>Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS FORMULATED STANDARDS SETS</b>	
<b>Subcategory: LEAD STANDARDS</b>	
<b>Table: IMPBT7A</b>	
<b>Standard No.</b>	<b>Lead (Grams/Gal.)</b>
1	0.000
2	0.100
3	1.000
4	2.000
5	3.000
6	4.000
7	5.000

<b>Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS FORMULATED STANDARDS SETS</b>	
<b>Subcategory: LEAD STANDARDS</b>	
<b>Table: IMPBT7C</b>	
<b>Standard No.</b>	<b>Lead (Grams/Gal.)</b>
1	0.000
2	0.001
3	0.005
4	0.010
5	0.050
6	0.100
7	0.300