

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: METAL WORKING FLUIDS**

**Table: METAL WORKING FLUIDS**

<b>Product Code</b>	<b>Product</b>	<b>Description</b>	<b>Price \$ Set</b>
IMMWFL13	METAL WORKING FLUIDS LUBRICANT Cl,P,S	For the analysis of Lubricant Base Metal Working Fluids. 13 standards per set, 100mls each	\$1,350.00

**Category: SULFUR, CHLORINE, LEAD OR LUBE OIL STANDARDS  
FORMULATED STANDARDS SETS**

**Subcategory: METAL WORKING FLUIDS**

**Table: IMMWFL13**

<b>Standard No.</b>	<b>Wt. % Cl</b>	<b>Wt. % P</b>	<b>Wt. % S</b>
1	0.0000	0.0000	0.0000
2	0.7500	0.0250	0.5000
3	0.0500	0.1000	3.0000
4	1.0000	0.5000	2.5000
5	0.1000	0.0050	2.0000
6	1.5000	0.2000	1.0000
7	2.0000	0.0050	3.0000
8	1.0000	0.0500	0.1000
9	0.5000	0.4000	0.0000
10	2.0000	0.2000	1.5000
11	0.0000	0.5000	1.5000
12	1.2500	0.0100	0.0500
13	0.0500	0.3000	0.0500