

Category: MONITOR SAMPLES FOR XRF - SILICATE GLASSES
Optically polished on the analytical surface; sizes 30-40mm in diameter and 5mm in height.

Subcategory: SINGLE ELEMENTS

Table: SINGLE ELEMENTS

Element	Product Code	Description	40 mm	Other Sizes
CHLORINE	IMGL-CL	Chlorine conc. 0.6wt%	\$275.00	\$295.00
SULFUR (HL)	IMGL-S(HL)	High level of Sulfur; approx. 0.70 wt%	\$275.00	\$295.00
SULFUR (LL)	IMGL-S(LL)	Low level of Sulfur; approx. 0.005 wt%	\$275.00	\$295.00
SULFUR (ML)	IMGL-S(ML)	Mid level of Sulfur; approx. 0.050 wt%	\$275.00	\$295.00

Category: MONITOR SAMPLES FOR XRF - SILICATE GLASSES
Optically polished on the analytical surface; sizes 30-40mm in diameter and 5mm in height.

Subcategory: MULTI-ELEMENTS

Table: Multi-Elements Price

Product Code	40 mm	Other Sizes
IMGL-LOE1	\$275.00	\$295.00
IMGL-LOE2	\$275.00	\$295.00
IMGL-LOE3	\$275.00	\$295.00

Category: MONITOR SAMPLES FOR XRF - SILICATE GLASSES
Optically polished on the analytical surface; sizes 30-40mm in diameter and 5mm in height.

Subcategory: MULTI-ELEMENTS

Table: Multi-Elements

Product Code Elements (approx wt%)	IMGL-LOE1	IMGL-LOE2	IMGL-LOE3
Ba	0.0040	0.1850	0.5600
Ca	0.0100	0.5000	0.5000
Cl	0.0050	0.1500	0.4500
Cu	0.0010	0.0500	0.1500
Mg	0.0010	0.1500	0.4500
Mo	0.0200	0.1500	0.4500
P	0.0000	0.1500	0.4600
S	0.2000	0.7500	0.7500
Zn	0.0200	0.1500	0.4600