

Category: STANDARDS FOR PVC ANALYSIS

Subcategory: PVC POWDERS

RoHS/WEEE compliant polymer standards in powder form for analysis of materials in polyvinylchloride (PVC). PVC's also include some Ca and Cl in each sample.

Available in 25 grams or 50 grams.

Table: PVC POWDERS

Product Code	Description	Price
IMPL(PVC)3-5E(P)	25 grams each and 1 QC Sample	\$650.00
IMPL(PVC)3-5E(P)	50 grams each and 1 QC Sample	\$1,000.00
IMPL(PVC)9-5E(P)	25 grams each and 1 QC Sample	\$1,750.00
IMPL(PVC)9-5E(P)	50 grams each and 1 QC Sample	\$2,650.00

Additional QC Samples

100 grams	\$400.00
250 grams	\$795.00
500 grams	\$1,100.00

Category: STANDARDS FOR PVC ANALYSIS

Subcategory: PL(PVC)3-5E(P)

Table: IMPL(PVC)3-5E(P)

Standard No.	Wt. % Br	Wt. % Cd	Wt.% Cr	Wt. % Hg	Wt. % Pb
1	0.0000	0.0000	0.0000	0.0000	0.0000
2	0.0500	0.0050	0.0500	0.0500	0.0500
3	0.1000	0.0100	0.1000	0.1000	0.1000

Wt. % Ca and Wt.% Cl will be 1.5000 and 15.00 in all 3 standards.

Category: STANDARDS FOR PVC ANALYSIS					
Subcategory: PL(PVC)9-5E(P)					
Table: IMPL(PVC)9-5E(P)					
Standard No.	Wt. % Br	Wt. % Cd	Wt.% Cr	Wt. % Hg	Wt. % Pb
1	0.0000	0.0000	0.0000	0.0000	0.0000
2	0.1000	0.0100	0.0010	0.0100	0.1000
3	0.1000	0.0100	0.1200	0.0150	0.0100
4	0.0100	0.0150	0.0750	0.0001	0.0010
5	0.0500	0.0075	0.1000	0.0075	0.1500
6	0.1500	0.0010	0.0500	0.0010	0.0075
7	0.0750	0.0001	0.0100	0.0200	0.0100
8	0.0100	0.0050	0.0500	0.0050	0.0500
9	0.0500	0.0150	0.0100	0.0050	0.0500
QC Sample	0.0500	0.0050	0.0500	0.0500	0.0500

Wt.% Ca will be as follows: 0, 0.500, 1.500, 0.250, 0, 3.0, 0, 2.0, 0 and 0 in the QC Sample.

Wt.% Cl will be as follows: 40.0, 35.0, 15.0, 40.0, 35.0, 10.0, 40.0, 12.5, 35.0 and 35.0 in the QC Sample.

Category: STANDARDS FOR PVC ANALYSIS

Subcategory: PVC DISC

RoHS/WEEE compliant polymer standards in disc form for analysis of materials in polyvinylchloride (PVC). Available in 40mm x 10mm (Chemplex Cups 1340) or 31mm x 10mm (Chemplex Cup 1330) disc with 1 Quality Control Sample. Supplied with the same concentrations as the PVC Powders.

Table: PVC DISC

Product Code	Description	Price
IMPL(PVC)3-5E(D)	31mm disc and 1 QC Sample	\$600.00
IMPL(PVC)3-5E(D)	40mm disc and 1 QC Sample	\$650.00
IMPL(PVC)9-5E(D)	31mm disc and 1 QC Sample	\$1,750.00
IMPL(PVC)9-5E(D)	40mm disc and 1 QC Sample	\$1,750.00

We also have a PE Additive Package that includes 5 standards with 10 elements each. The elements included are Al, Ca, F, Mg, Na, P, S, Si, Ti and Zn.