

Instrumed

Instrumed Laboratory

Chemistry Analyzer



Abbott Istat

[Call for Pricing](#)

The i-STAT System represents the next step in the evolution of physician office testing. It allows your office to minimize treatment and follow-up delays that can result from sending patients to an off-site lab testing facility. The i-STAT System deliv [\[Product Details...\]](#)



ATAC 8000

[Call for Pricing](#)

The ATAC 8000, is a Random Access Chemistry System for all of your Clinical Chemistry Testing. [\[Product Details...\]](#)



Ortho Vitros 250/350

[Call for Pricing](#)

The VITROS® 250 and 350 Chemistry System offers easy operation, maintenance, and training. Whether you use it as a primary, STAT, or back-up to the powerful VITROS® 5,1 FS, the VITROS® 250 and 350 Chemistry System will get the job done right the first [\[Product Details...\]](#)



Roche Cobas Mira Plus CC

[Call for Pricing](#)

The Roche Cobas Mira is a benchtop, random access analyzer. This system allows for the selective analysis of chemistries in either a routine or STAT mode. Testing sequence is by patient rather than test (batch). The optional ISE module for sodium, potassi [\[Product Details...\]](#)



Roche Hitachi 902 Chemistry Analyzer

[Call for Pricing](#)

Reliability, quality and convenience make the Roche/Hitachi 902 the clinical chemistry workhorse for the small workload laboratory. Compact, sample selective clinical chemistry analyzer with a capacity up to 36 tests on board (including 3 ISE). [\[Product Details...\]](#)



Roche Hitachi 911 Chemistry Analyzer

[Call for Pricing](#)

A fully automated, discrete analyzer comprising of two units: control unit and analytical unit. Control unit is composed of a keyboard, a monitor (CRT), and a printer. The analytical unit involves a photometric measuring system, an ISE System and a CPU. [\[Product Details...\]](#)



Roche Hitachi 912 Chemistry Analyzer

[Call for Pricing](#)

Reliability, quality and convenience make the Roche/Hitachi 912 the optimal clinical chemistry workhorse system for the small up to medium workload laboratory. Sample selective analyzer designed for mid-size workloads. [\[Product Details...\]](#)



Roche Hitachi 917 Chemistry Analyzer

[Call for Pricing](#)

The Hitachi 917 analyzer is a fully automated, discrete, computerized chemistry analyzer. It is intended for in vitro quantitative or qualitative determination of a wide range of analytes in blood and other body fluids.

[\[Product Details...\]](#)



Roche Integra 400+

[Call for Pricing](#)

Integra 400+: The pressure for laboratories to work more efficiently has never been higher, increasing the need for innovative solutions. With the COBAS INTEGRA® 400 plus, innovation means 4 on-board measurement technologies, robotic handling of reaction [\[Product Details...\]](#)



Roche Integra 800

[Call for Pricing](#)

Integra 800: Sample selective bench top analyzer for small to mid-size laboratories. Flexible system for consolidating routine workload with a capacity of 36 tests (including 4 ISE: Na+, K+, Li+, Cl-) on board and a throughput up to 400 test per hour. [\[Product Details...\]](#)



Siemens Advia 1650 Chemistry Analyzer

[Call for Pricing](#)

The ADVIA® 1650 has minimal requirements for sample volume (as low as 2 μ L/test), reagent consumption (as low as 80 μ L/test) and water consumption (25 L/hour) - thus maximizing management efficiency. A rapid, 3-second cycle and system throughput of up [\[Product Details...\]](#)



Siemens Dimension RXL MAX

[Call for Pricing](#)

* Ability to run up to 91 methods onboard concurrently * True integration of chemistry and immunoassay with uncompromised throughput * Built-in automation provides the confidence that anyone can run any test at any time * Full disease-state [\[Product Details...\]](#)



Siemens Dimension Xpand Plus

[Call for Pricing](#)

* Full-range integrated capability in the smallest footprint for an instrument of its capabilities * True integration of chemistry and immunoassay provides uncompromised throughput * Built-in automation provides the confidence that anyone can run [\[Product Details...\]](#)



Vitros DT60II Chemistry Analyzer

[Call for Pricing](#)

Get all the benefits of VITROS® MicroSlide® Technology in a size that's perfect for your laboratory. Fast, more accurate results with an easy-to-use, cost-effective solution. The modular VITROS® DT60 II Chemistry System suits your workflow perfectly. See [\[Product Details...\]](#)

- « « Start
- « Prev
- 1
- Next »
- End » »

Results 1 - 14 of 14