

# Instrumed

*Instru-Med Scientific*

## Consumables



### Disposable Sample Cells

[Call for Pricing](#)

To ensure extremely easy, reliable analysis, HORIBA has developed a unique, disposable sample cell, allowing you to easily see the sample in a leak-free package without spilling any of the sample. [\[Product Details...\]](#)



### SOLVENT S-316

[Call for Pricing](#)

The environmentally friendly S-316 extraction solvent is used with the OCMA-350 to extract oil components from oily water samples, soil samples, or product surfaces. [\[Product Details...\]](#)

- - 
  - 
  - 
  -
- « « Start  
« Prev  
1  
Next »  
End » »

Results 1 - 2 of 2