

Instrumed

Instru-Med Scientific

AGILITY PORTABLE BENCHTOP RAMAN SPECTROMETER

[Call for Pricing](#)



BaySpec's Agility Raman spectrometer delivers high sensitivity and repeatability in an affordable, ruggedized, battery-operated package. Available in 532, 785, and 1064 nm with single or dual band options. An integrated sample compartment allows the utmost flexibility, via its quick-change, auto-aligning sample holders.

Features:

- Integrated laser, spectrometer, touch screen computer
- Standard wavelength options at 532, 785, and 1064 nm (custom wavelengths available)
- Dual wavelength options at 532/1064 and 785/1064 nm (custom wavelengths available)
- Various Quick-Change Sample options for any sample type
- Cooled CCD or InGaAs detector
- Portable ergonomic design, light-weight
- High throughput transmission spectrograph with VPG and f/1.8 design
- Up to 499 mW laser power at 785 and 1064 nm, up to 50 mW at 532 nm
- Agile20/20 Windows-based software for turn-key operation

For more product information please visit:

<http://www.bayspec.com/spectroscopy/agility-portable-benchtop/>

[Vendor Information](#)

